

## Installing CancerHelp™ Monthly CD

### Quit from the CancerHelp™ Program

First connect the keyboard and mouse to the blue end of the “Y” splitter or to the keyboard extension box, if it is not already attached. If CancerHelp™ is running, click the mouse button or touch the screen once to get out of the screen saver and back to the “English or Spanish” Screen.



While holding down the Apple key on the keyboard (just to the left of the space bar), press and release the letter “Q” on the keyboard. This will Quit out of CancerHelp™.

You should see the Desktop Screen like that shown below:



### Insert the CancerHelp™ CD-Rom

On the front of the CPU is a button that opens a CD-Rom tray. Open the tray and insert the CancerHelp™ monthly CD. After a few seconds the icon for the CD-Rom will appear in the upper right side of the screen just below the icon for “Macintosh HD” or “Hard Disk.”



After you see the icon for the CD-Rom, you should also see the window for the CD-Rom in the center of the screen.



You will see an icon with a hand on it and text “CancerHelp dbl click.” Double click on this icon to launch the CancerHelp™ Monthly Installer.

## Installing CancerHelp™ Monthly CD

### Running the CancerHelp™ Installer

After you double click on the icon, a window will open with the text “CancerHelp™ Monthly Installer” and a button with the word “Install.”



Click once on the “Install” button and the Installer will start to copy the files for this month’s version to the hard drive of your computer. Specifically they are being copied to a folder called “CancerHelp™ Program.”

When the Installer finishes, a window will appear with “Congratulations...” and the following information.



Click the “Quit” button in the lower right corner of this text window.

### Closing the Window & Ejecting the CD

Now the screen looks like it did when you first put in the CD-Rom.



Close the CD-Rom Window (the one with the “CancerHelp™ Install.dbl click” icon in it) by taking the mouse arrow and clicking on the small box that is in the upper left corner of the CD-Rom Window.


Ejecting the CD-Rom from the computer is a little tricky so you may need to try it several times. Take the mouse arrow to on top of the CD-Rom icon in the upper right corner of the screen just below the icon for “Hard Disk” or “Macintosh HD.”

Click once on the icon and while holding the mouse button down slide the icon for the CD-Rom to on top of the trash can icon until you see the trash can highlight or turn black. When you see the trash can highlight, let go of the mouse and the CD-Rom icon will disappear and the CD-Rom tray will open. Take out the CD-Rom and close the tray. Keep two months of CD-Roms to reinstall CancerHelp™ if needed.

## Installing CancerHelp™ Monthly CD

### Restarting CancerHelp™

When the CD-Rom has ejected, you can restart CancerHelp™ by restarting the computer. To Restart the computer, take the mouse pointer to the “Special” menu at the top of the screen. (Note the menu bar has an Apple icon and then the words File, Edit, View and Special in the top left corner of the screen. After selecting restart the computer will



relaunch and after a minute or so the CancerHelp program will start up.

### Testing the Touchscreen Monitor and Printer

After CancerHelp reaches the English or Spanish screen, test the touchscreen by navigating through the software. Observe whether where you touch on the screen is where it responds. If the response is off, then you can recalibrate the touchscreen monitor. Also find a topic to print and print to the printer. Check that there is enough paper in the paper tray and the ink on the print out is ledgeble and clear to read. If the print out is faint, replace the toner or ink cartridge.